


## Cleaning brushes


### Cleaning brushes



Autoclavable

Code	Article	Description	Length (mm)	Outer Diameter (mm)
82287	CS-3025S	Cleaning brushes for L/W - and WJ-Channel	2500	3.0
82479	CS-3030S	Cleaning brushes extra long, for L/W - and WJ-Channel	3000	3.0
82288	CS-3010S	Cleaning brushes working channel with max. 1.2 mm diameter	1000	3.0
82292	CS-6002SN	Cleaning brushes working channel, short	200	6.0
81454	CS-6015ST	Cleaning brushes for working channel	1500	6.0
82290	CS-6021SN	Cleaning brushes for working channel	2100	6.0
82289	CS-6030SN	Cleaning brushes for working channel, extra long	3000	6.0

### Cleaning brushes for single use only



For single use only

Code	Article	Description	Length (mm)	Outer Diameter (mm)
86151	CS-3045	Single-use cleaning brushes for WJ channel, for 70(K), 80K, 81K, 10 pcs.	4500	3.0
86147	CS-6021T	Single-use cleaning brushes for working channel for 70(K), 80K, 81K, 85(K), 90K, 90i series, 10 pcs.	2100	6.0

**PENTAX Europe GmbH**  
LIFE CARE  
Julius-Vosseler-Straße 104  
22527 Hamburg  
Germany  
Tel.: +49 40 / 5 61 92 - 0  
Fax: +49 40 / 5 60 42 13  
E-Mail: medical@pentax.de

**PENTAX U.K. Limited**  
LIFE CARE  
Pentax House  
Heron Drive, Langley  
Slough SL3 8PN  
United Kingdom  
Tel.: +44 17 53 / 79 27 23  
Fax: +44 17 53 / 79 27 94  
E-Mail: lifecare.info@pentax.co.uk

**PENTAX France S.A.S.**  
LIFE CARE  
112 quai de Bezons  
P. B. 204  
95106 Argenteuil  
France  
Tel.: +33 1 / 30 25 75 75  
Fax: +33 1 / 30 25 75 76  
E-Mail: medical@pentax.fr

**PENTAX Nederland B.V.**  
LIFE CARE  
Lage Mosten 35  
4822 NK Breda  
Netherlands  
Tel.: +31 76 / 5 31 30 31  
Fax: +31 76 / 5 31 30 00  
E-Mail: medical@pentax.nl

**PENTAX Italia S.r.l.**  
LIFE CARE  
Via Dione Cassio, 15  
20138 Milano  
Italy  
Tel.: +39 / 02 50 99 58 1  
Fax: +39 / 02 50 99 58 60  
E-Mail: marketing.lifecare@pentaxitalia.it

**S.I.M. - Sitemas Integrales de Medicina, S.A.**  
LIFE CARE  
Calle Mario Cabré, 13  
28030 Madrid  
Spain  
Tel.: +034-91-3016240  
Fax: +034-91-7513115

**HOYA Corporation**  
**PENTAX Tokyo office**  
2-36-9, Maeno-cho  
Itabashi-Ku  
174-8639 Tokyo  
Japan  
Tel.: +81 33 / 9 60 51 55  
Fax: +81 35 / 3 92 67 24

# PENTAX

## Szczoteczki czyszczące

### Szczoteczki czyszczące

/rysunek/ /rysunek/ autoklawowalne

kod	artykuł	opis	Długość [mm]	Średnica zewnętrzna [mm]
82287	CS-3025S	Szczoteczka czyszcząca do kanałów L/W i kanału WJ	2500	3,0
82479	CS-3030S	Szczoteczka czyszcząca bardzo długa, odpowiednia do kanałów L/W i kanału WJ	3000	3,0
82288	CS-3010S	Szczoteczka czyszcząca do kanału roboczego o max. średnicy 1,2 mm	1000	3,0
82292	CS-6002SN	Szczoteczka czyszcząca, krótka	200	6,0
81454	CS-6015ST	Szczoteczka czyszcząca do kanału roboczego	1500	6,0
82290	CS-6021SN	Szczoteczka czyszcząca do kanału roboczego	2100	6,0
82289	CS-6030SN	Szczoteczka czyszcząca do kanału roboczego, bardzo długa	3000	6,0

### Szczoteczki czyszczące jednorazowego użytku

/rysunek/ /rysunek/ jednorazowego użytku

kod	artykuł	opis	Długość [mm]	Średnica zewnętrzna [mm]
86151	CS-3045	Szczoteczka czyszcząca jednorazowego użytku do kanału WJ, do 70(K), 80K, 81K, 10 szt	4500	3,0
86147	CS-6021T	Szczoteczka czyszcząca jednorazowego użytku do kanału roboczego do 70(K), 80K, 81K, 85(K), 90K, serii 90i, 10 szt	2100	6,0

**PENTAX Europe GmbH**  
LIFE CARE  
Julius-Vosseler-Straße 104  
22527 Hamburg  
Germany  
Tel.: +49 40 / 5 61 92 - 0  
Fax: +49 40 / 5 60 42 13  
E-Mail: medical@pentax.de

**PENTAX U.K. Limited**  
LIFE CARE  
Pentax House  
Heron Drive, Langley  
Slough SL3 8PN  
United Kingdom  
Tel.: +44 17 53 / 79 27 23  
Fax: +44 17 53 / 79 27 94  
E-Mail: lifecare.info@pentax.co.uk

**PENTAX France S.A.S.**  
LIFE CARE  
112 quai de Bezons  
P. B. 204  
95106 Argenteuil  
France  
Tel.: +33 1 / 30 25 75 75  
Fax: +33 1 / 30 25 75 76  
E-Mail: medical@pentax.fr

**PENTAX Nederland B.V.**  
LIFE CARE  
Lage Mosten 35  
4822 NK Breda  
Netherlands  
Tel.: +31 76 / 5 31 30 31  
Fax: +31 76 / 5 31 30 00  
E-Mail: medical@pentax.nl

**PENTAX Italia S.r.l.**  
LIFE CARE  
Via Dione Cassio, 15  
20138 Milano  
Italy  
Tel.: +39 / 02 50 99 58 1  
Fax: +39 / 02 50 99 58 60  
E-Mail: marketing.lifecare@pentaxitalia.it

**S.I.M. - Sitemas Integrales de Medicina, S.A.**  
LIFE CARE  
Calle Mario Cabré, 13  
28030 Madrid  
Spain  
Tel.: +034-91-3016240  
Fax: +034-91-7513115

**HOYA Corporation**  
**PENTAX Tokyo office**  
2-36-9, Maeno-cho  
Itabashi-Ku  
174-0639 Tokyo  
Japan  
Tel.: +81 33 / 9 60 51 55  
Fax: +81 35 / 3 92 67 24

# PENTAX